AED Monthly Maintenance Checklist

Refer to your AED manufacturer's maintenance recommendations

School/Location:							SY:	20					
Defibrillator Unit Make/Model:					Unit Serial Number:								
*AED Unit	Batteries are replaced every (<i>X)</i>	years depending on your AED manufacturer's recommendations											
Instruction	Recommended / Corrective Action	Sept	Oct	Nov	Dec	Jan	Feb	Mar	Apr	May	June	Jul	Aug
1. Check readiness display for: *Insert specific characteristics of school AED unit * Example: Green check mark indicator 2. Check and note expiration date on electrode pads: *Adult: *Peds (if applicable): 3. Check additional supplies: *Extra Adult Electrodes *Pouch: scissors, blob (mask), razor, gloves, antiseptic towelette	Replenish as needed												
4. Examine the AED, case, cables and connectors for: *Damage or Cracks *Foreign Substances	*Example: Contact authorized service provider Clean the device with soap and damp cloth												
5. Wall cabinet: *Case is intact, no tampering has been noted *Alarm battery (if applicable)	*Example: 9-volt battery to be changed annually in June, prior to the end of school year by custodial staff												
AED Maintenance Inspectors: Name / Initials			Name /	Initials									